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Visit our website: www.dentsplymaillefer.com

FOR DENTAL USE ONLY

GB DIRECTIONS FOR USE (Thermanprep[®] 2 oven)

- 1) Indications for use**
This product has to be used only in a clinical or hospital environment, by qualified users.
The Thermanprep[®] 2 Obturator Oven has been specially developed for heating Dentsply Maillefer endodontic obturators (GuttaCore™, WaveOne™, Thermanfil[®], Protaper™, GT Series X™) for use in root canal therapy. Do not use this oven for any other purpose.
- 2) Contraindications**
There are none known.
- 3) Warnings**
To abort the cleaning mode without disconnecting the main plug and to enter the normal mode again press either button for 2 seconds.
- 4) Precautions**
 - To prevent fire or electric shock, do not expose this product to moisture.
 - Unplug the oven when not in use.
 - Do not remove the cover of the oven.
 - This oven is designed for use in the dental office. Do not use in any other environment.
 - To reduce the risk of electric shock, do not remove cover, there are no user-serviceable parts inside.
 - This oven is intended for use in heating gutta-percha coated obturators only. Do not insert anything other than an obturator into the oven.
 - Changes or modifications not expressly approved by DENTSPLY Tulsa Dental Specialties could void the user's authority to operate the equipment.
- 5) Adverse reactions**
In the present technical state, no adverse reaction has been reported so far.
- 6) Step-by-step instructions**
 - Turn the power on by plugging the oven into a 115 or 230 VAC outlet. Both indicator lights around the buttons will light up. To wake-up the oven from stand-by mode, push and hold either one of the buttons until the indicator light around the button is illuminated and two signals "beeps" are heard, and release the button. A second, lower signal "beep" will sound. On both sides the smallest obturator size is selected as shown by indicator lights. The oven is now ready for use. At any time the oven can be turned off. Push and hold either one of the bottom buttons until the indicator light around the button turns off.
 - Put both obturator holders in the upper position by pressing the back of the holder (marked area).
 - DISINFECT the obturator in a 5.25% sodium hypochlorite solution (bleach) for one minute, RINSE the obturator in 70% alcohol. DRY the obturator for several seconds on a clean work surface to allow evaporation of excess alcohol.
 - After disinfecting the obturator and setting the rubber stop, place the Obturator in one of the obturator holders, lightly press on top to have it firmly seated. Regardless of what working length you have determined, the rubber stop must be under the holder. Make sure the obturator is hanging straight, not at an angle. The obturator must clear all sides of the heating chamber.
 - Push the button which corresponds to the selected obturator holder, to select the size of the Obturator you wish to heat. Repeatedly push this button to select the correct size.
 - To start the heating process push the obturator holder down until you hear it click. It will raise a little but keep in a steady down position. The indicator light around the button will show that the heater is on (constantly on). The obturator is now being heated. Heating time varies – depending on obturator size and type – from 20 to 49 seconds, and is regulated automatically.
The heating times for the selected obturator sizes are as follows:

Use the following guide for heating Dentsply Maillefer endodontic obturators:

OBTURATOR BRAND	STICKER	OBTURATOR SIZE	HEATING POSITION	SIGNAL TIME	TOLERANCE
GUTTA CORE	GUTTA CORE	20 to 90	20 - 90	20 seconds	± 5 seconds
		small	Small / Primary	20 seconds	± 5 seconds
		primary	primary	20 seconds	± 5 seconds
		large	large	44 seconds	± 7 seconds
		F1	F1-F2	20 seconds	± 5 seconds
WAVE ONE	WAVE ONE	F2	F2	20 seconds	± 5 seconds
		F3	F3	44 seconds	± 7 seconds
		F4	F3-F4-F5	44 seconds	± 7 seconds
		F5	F5	20 seconds	± 5 seconds
		20-25	20-25	20 seconds	± 5 seconds
THERMAFIL	THERMAFIL	30-35-40-45	30-60	44 seconds	± 7 seconds
		50-55-60	30-60	44 seconds	± 7 seconds
		70-80-90-100	70-140	49 seconds	± 8 seconds
GT series X	GT series X	110-120-130-140	70-140	49 seconds	± 8 seconds
		.04	Size .04-.08	20 seconds	± 5 seconds
		.06	Size .04-.08	20 seconds	± 5 seconds
		.08	Size .04-.08	20 seconds	± 5 seconds
		.12	Size .10-.12	44 seconds	± 7 seconds

- After the first signal "beep", the obturator is ready for use. Push the back of the obturator holder (marked area) to let it raise slowly. The indicator light around the button will go off. Take the obturator carefully out of the holder by pulling it towards you, making sure not to scrape the obturator on any part of the holder.
- You may leave the obturator in the Thermanprep[®] 2 Obturator Oven for up to 90 seconds after the first signal "beep". It will keep the obturator at the appropriate temperature and ready for use. The oven will "beep" every 15 seconds to remind you that the obturator is still in the oven. After 90 seconds, the heating element will switch off automatically.
- If you want to heat more obturators, you may use one holder or both to continue your work efficiently. Each holder is operated by its own button and indicators. If you hear a signal "beep" then look at indicator lights around the buttons which indicate obturator heater is ready (flashing). When both heaters have finished their cycle and no buttons are pushed, the oven will go into a power-saving, stand-by mode after 5 minutes.
- See below table for oven specifications:

The equipment is designed to be safe for indoor use under following environmental conditions:	
Max. altitude	2 000 m
Storage temperature	-20 to 65°C -4 to +149°F
Storage humidity	20 to 80% no condensation
Operating temperature	+5 to +40°C +41 to +104°F
Max. operating humidity	80% at +31°C to 50% at 40°C (lin.) 80% at +87°F to 50% at 122°F (lin.)
Mains voltage fluctuations	Up to ±10% of the nominal voltage
Transient overvoltages	Category II
Rated pollution degree	2
Mass of the equipment	0.56 kg
Overall size of the equipment (W x D x H)	155 x 178 x 81 mm

If there is any doubt about the compatibility of decontamination or cleaning agents, please contact the cleaning agent's manufacturer or DENTSPLY Tulsa Dental Specialties

NOTE: Precision Heating Alerts

The oven is designed with two precision heating alerts:

- First size LED button blinks:** The targeted temperature of the heating element was not reached within specified time (time-out). Use the other heating element to warm the obturator. This may be re-set by turning the power off and turning the power on. If the error occurs again the oven should no longer be used and sent for repair.
- Second size LED button blinks:** The hardware safety circuit has been activated due to a failure in the heating element, causing overheating. This may be re-set by turning power off, disconnect the power cord, reconnect the power cord, and turn the power on. If the error occurs again the oven should no longer be used and sent for repair.

7) Cleaning instructions

Clean oven after each use. Any gutta-percha remaining on the holder or in the heating chamber may be effectively removed with Orange Solvent (generic dental solvent) on a cotton pledget. Press and hold both buttons at once for 2 seconds to put the oven into the cleaning mode. This will cause both heaters to heat up to a medium temperature of 140°F/60°C. The indicator lights around the buttons will show that the heaters are on (constantly on). A signal "beep" and flashing indicator lights around the buttons will warn that the desired temperature has been reached. Disconnect the main plug and clean the heater chambers with a wire brush from bottom side as indicated.

8) Technical specifications


Voltage rating: 115 or 230 VAC, 50/60 Hertz, range of the supply voltage ±10%, power-save stand-by mode: <1.5W, normal use: 380W max Fuse Rating: 4AT 250V

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Printed on Thu Jun 16 2016 17:41:26 GMT+0300

Explanation of Marking Symbols	
	Recycling: Please do not throw away! This product and all its components must absolutely be recycled through your supplier.
	Manufacturer
	See directions for use
	Consult instructions for use
	Opened packages are not replaced
	Class II equipment
	Hot surface
	Serial Number

- After selecting proper size press here to lower the obturator holder and start heating Press here again to stop heating and raise the obturator holder
- Place obturator here
- Heating chamber
- Sizes 70-140 / 30-60 / 20-25
- Press button to select obturator size for left obturator heater
- LEFT heater "ON" indicator light
- Cleaning mode indicator light
- Press button to select obturator size for right obturator heater
- RIGHT heater "ON" indicator light
- Main cord, disconnect after use
- Clean heating chambers from here with supplied tool
- WARNING: Disconnect main power cord before entering heating chambers

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Dentsply De Trey GmbH